# antibodies -online.com





#### Datasheet for ABIN7443398

## anti-HVD3 antibody



#### Overview

Quantity:	100 μL
Target:	HVD3
Reactivity:	Various Species
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This HVD3 antibody is un-conjugated
Application:	Western Blotting (WB), Immunohistochemistry (IHC), Immunoprecipitation (IP), Immunocytochemistry (ICC)

#### **Product Details**

Target:

Purpose:	Polyclonal Antibody to 25-Hydroxyvitamin D3 (HVD3)
'	
Immunogen:	OVA Conjugated 25-Hydroxyvitamin D3 (HVD3)
Isotype:	IgG
Specificity:	The antibody is a rabbit polyclonal antibody raised against HVD3. It has been selected for its
	ability to recognize HVD3 in immunohistochemical staining and western blotting.
Cross-Reactivity:	Various Species
Purification:	Antigen-specific affinity chromatography followed by Protein A affinity chromatography
Target Details	

HVD3

#### Target Details

Alternative Name:	25-Hydroxyvitamin D3 (HVD3 Products)
Target Type:	Chemical
Background:	25-HVD3, 25HVD3, Calcifediol, Calcidiol, H-VD3, 25-Hydroxycholecalciferol, 25-Hydroxy Vitamin D3, 25(OH)D

### **Application Details**

Application Notes:	Immunohistochemistry: 5-20 μg/mL
	Immunocytochemistry: 5-20 μg/mL
	Optimal working dilutions must be determined by end user.
Comment:	The thermal stability is described by the loss rate. The loss rate was determined by accelerated
	thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious
	degradation and precipitation were observed. The loss rate is less than 5% within the expiration
	date under appropriate storage condition.
Restrictions:	For Research Use only

#### Handling

Format:	Liquid
Buffer:	0.01M PBS, pH 7.4, containing 0.05 % Proclin-300, 50 % glycerol.
Preservative:	ProClin
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only.
Storage:	4 °C,-20 °C
Storage Comment:	Store at 4°C for frequent use. Stored at -20°C in a manual defrost freezer for two year without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Expiry Date:	24 months