
antibodies .- online.com







anti-PIIINP antibody

Images



Overview

Quantity:	100 μL
Target:	PIIINP
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This PIIINP antibody is un-conjugated
Application:	Western Blotting (WB), Immunohistochemistry (IHC), Immunoprecipitation (IP), Immunocytochemistry (ICC)

Product Details


Alternative Name:

1 Toddet Details	
Purpose:	Polyclonal Antibody to Procollagen III N-Terminal Propeptide (PIIINP)
Immunogen:	OVA Conjugated Procollagen III N-Terminal Propeptide (PIIINP)
Isotype:	IgG
Specificity:	The antibody is a rabbit polyclonal antibody raised against PIIINP. It has been selected for its ability to recognize PIIINP in immunohistochemical staining and western blotting.
Purification:	Antigen-specific affinity chromatography followed by Protein A affinity chromatography
Target Details	
Target:	PIIINP

Procollagen III N-Terminal Propeptide (PIIINP Products)


Target Details

3 - 1 - 1	
Background:	P3NP, N-Propeptide Of Type III Procollagen, Procollagen III Amino Terminal Propeptide
Application Details	
Application Notes:	Western blotting: 0.5-2 μg/mL
	Immunohistochemistry: 5-20 μg/mL
	Immunocytochemistry: 5-20 μg/mL
	Optimal working dilutions must be determined by end user.
Comment:	The thermal stability is described by the loss rate. The loss rate was determined by accelerated
	thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious
	degradation and precipitation were observed. The loss rate is less than 5% within the expiration
	date under appropriate storage condition.
Restrictions:	For Research Use only
Handling	
Format:	Liquid
Buffer:	PBS, pH 7.4, containing 0.02 % Sodium azide, 50 % glycerol.
Preservative:	Sodium azide
Precaution of Use:	This product contains Sodium azide: a POISONOUS AND HAZARDOUS SUBSTANCE which
	should be handled by trained staff only.
Storage:	4 °C,-20 °C
Storage Comment:	Store at 4°C for frequent use. Stored at -20°C in a manual defrost freezer for two year without
	detectable loss of activity. Avoid repeated freeze-thaw cycles.
Expiry Date:	24 months



Immunohistochemistry

Image 1. Detection of PIIINP in Human Pancreas Tissue using Polyclonal Antibody to Procollagen III N-Terminal Propeptide (PIIINP)



Immunohistochemistry

Image 2. Detection of PIIINP in Human Cardiac Muscle Tissue using Polyclonal Antibody to Procollagen III N-Terminal Propeptide (PIIINP)