



[Go to Product page](#)

Datasheet for ABIN7453138  
**anti-NASP antibody (AA 150-200)**

### Overview

|                      |                                     |
|----------------------|-------------------------------------|
| Quantity:            | 20 µg                               |
| Target:              | NASP                                |
| Binding Specificity: | AA 150-200                          |
| Reactivity:          | Human                               |
| Host:                | Rabbit                              |
| Clonality:           | Polyclonal                          |
| Conjugate:           | This NASP antibody is un-conjugated |
| Application:         | Western Blotting (WB)               |

### Product Details

|               |  |
|---------------|--|
| Purpose:      | Rabbit anti-NASP Antibody, Affinity Purified |
| Immunogen:    | between AA 150 and 200                       |
| Isotype:      | IgG  |
| Purification: | Affinity Purified                            |

### Target Details

|                   |  |
|-------------------|--|
| Target:           | NASP   |
| Alternative Name: | NASP ( <a href="#">NASP Products</a> )   |
| Background:       | Background: NASP (Nuclear Autoantigenic Sperm Protein) is a component of a multichaperone nucleosome-remodeling complex that shuttles histones and coordinates their deposition in |

## Target Details

---

new nucleosomes. NASP functions as an H1 histone chaperone that transports H1 histones into the nucleus and exchanges H1 histones with DNA. NASP has been demonstrated to be required for DNA replication and cell cycle progression.

Gene ID: 4678

UniProt: [P49321](#)

## Application Details

---

Application Notes: IP: Not recommended  
WB: 1:2,000 - 1:10,000

Restrictions: For Research Use only

## Handling

---

Concentration: 200 µg/mL

Buffer: Tris-buffered Saline containing 0.1 % BSA and 0.09 % Sodium Azide

Preservative: Sodium azide

Precaution of Use: This product contains Sodium azide: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only.

Storage: 4 °C

Expiry Date: 12 months