

## Datasheet for ABIN7477525

# Recombinant anti-MUC1 antibody





$\sim$			
( )	ve	r\/	٨
$\cup$	V	1 V I	/١


Quantity:	1 mL	
Target:	MUC1	
Reactivity:	Human	
Host:	Rabbit	
Antibody Type:	Recombinant Antibody	
Clonality:	Monoclonal	
Conjugate:	This MUC1 antibody is un-conjugated	
Application:	Immunohistochemistry (IHC)	
Product Details		
Immunogen:	Recombinant human MUC1 protein	
Clone:	MSVA-672R	
Isotype:	IgG	
Target Details		
Target:	MUC1	
Alternative Name:	MUC1 / EMA (MUC1 Products)	

Mucin-1 subunit beta, Peanut-reactive urinary mucin, PEM, PEMT, Polymorphic epithelial mucin,

### **Target Details**


rarget Details		
	PUM, Tumor-associated epithelial membrane antigen,MUC1 / EMA antibody validated for	
	immunohistochemistry on 76 different Normal Tissues	
UniProt:	P15941	
Pathways:	Negative Regulation of intrinsic apoptotic Signaling	
Application Details		
Application Notes:	IHC 1:75-1:150	
Comment:	Positive Control: Kidney: Distal tubuli and collecting ducts must show a strong cytoplasmic MUC1 staining while the staining should be at least moderate in the parietal and visceral sheets of the Bowman's capsule and largely absent in proximal tubuli.  Negative Control: Colon: MUC1 staining should be very weak or absent in colonic epithelium.	
Protocol:	Manual Protocol: Freshly cut sections should be used (less than 10 days between cutting and staining). Heat-induced antigen retrieval for 5 minutes in an autoclave at 121 °C in pH 7,8 Target Retrieval Solution buffer. Apply the antibody at a dilution of 1:100 at 37 °C for 60 minutes. Visualization of bound antibody by the EnVision Kit (Dako, Agilent) according to the manufacturer's directions.	
Restrictions:	For Research Use only	
Handling		
Format:	Liquid	
Storage:	4 °C,-20 °C,-80 °C	
Storage Comment:	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody	

is stable for 24 months. Non-hazardous.




#### **Immunohistochemistry**

**Image 1.** Pancreas Pancreatic acinar cells show a strong MUC1 immunostaining of the apical membranes and of granular Golgi like structures below the apical surface



#### **Immunohistochemistry**

Image 2. Kidney cortex



#### **Immunohistochemistry**

**Image 3.** Malignant mesothelioma epitheloid variant with strong predominantly membranous EMA immunostaining of tumor cells