

Datasheet for ABIN749222

anti-Plexin A1 antibody (AA 901-1000) (HRP)



Overview

Overview	
Quantity:	100 μL
Target:	Plexin A1 (PLXNA1)
Binding Specificity:	AA 901-1000
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This Plexin A1 antibody is conjugated to HRP
Application:	ELISA, Immunohistochemistry (Frozen Sections) (IHC (fro)), Immunohistochemistry (Paraffinembedded Sections) (IHC (p))
Product Details	
Immunogen:	KLH conjugated synthetic peptide derived from human Plexin A1
Isotype:	IgG
Cross-Reactivity:	Mouse
Predicted Reactivity:	Human,Rat,Dog,Pig,Horse,Rabbit,Guinea Pig
Purification:	Purified by Protein A.
Target Details	
Target:	Plexin A1 (PLXNA1)
Alternative Name:	Plexin A1 (PLXNA1 Products)

Target Details

Background:	Synonyms: NOV, NOVP, PLXN1, PLEXIN-A1, Semaphorin receptor NOV, PLXNA1
	Background: Coreceptor for SEMA3A, SEMA3C, SEMA3F and SEMA6D. Necessary for signaling
	by class 3 semaphorins and subsequent remodeling of the cytoskeleton. Plays a role in axon
	guidance, invasive growth and cell migration. Class 3 semaphorins bind to a complex
	composed of a neuropilin and a plexin. The plexin modulates the affinity of the complex for
	specific semaphorins, and its cytoplasmic domain is required for the activation of down-stream
	signaling events in the cytoplasm (By similarity).
Gene ID:	5361
UniProt:	Q9UIW2
Pathways:	Smooth Muscle Cell Migration
Application Details	
Application Notes:	IHC-P 1:200-400
	IHC-F 1:100-500
Restrictions:	For Research Use only
Handling	
Format:	Liquid
Concentration:	1 μg/μL
Buffer:	Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.
Preservative:	ProClin
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.
Handling Advice:	Do NOT add Sodium Azide! Use of Sodium Azide will inhibit enzyme activity of horseradish peroxidase.
Storage:	-20 °C
Storage Comment:	Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.
Expiry Date:	12 months