

## Datasheet for ABIN7526137 Cell Adhesion Molecule 4 IQ-ELISA Kit



Overview

Quantity:	
Quantity.	2 x 96 tests
Target:	Cell Adhesion Molecule 4 (CADM4)
Reactivity:	Human
Method Type:	Sandwich ELISA
Application:	ELISA
Product Details	
Purpose:	Human Immunoquantitative (PCR-Based) IGSF4C ELISA Kit for cell culture supernatants,
	plasma, and serum samples.
Sample Type:	Cell Culture Supernatant, Plasma, Serum
Analytical Method:	Semi-Quantitative
Detection Method:	qPCR
Characteristics:	IQ-ELISAs employ specific capture antibodies coated on a 96-well PCR plate. Standards and
	samples are pipetted into the wells, the target protein in the standards and samples binds to
	the immobilized antibody. The wells are washed and the detection affinity reagent is added to
	the wells where it binds to any captured antigen. The wells are washed, and primers and PCR
	master mix are added to each well. The plate is placed into a real time PCR instrument for
	cycling and measurement of DNA amplification. The cycle number where amplification is
	detected is proportional to the amount of affinity detection reagent that bound to captured
	antigen in each well.
Components:	IGSF4C Microplate (Item A): 96 well PCR plate coated with anti-Human IGSF4C

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/3 | Product datasheet for ABIN7526137 | 01/06/2025 | Copyright antibodies-online. All rights reserved.

	Wash Buffer I Concentrate (20x) (Item B): 25 ml of 20x concentrated solution		
	Standards (Item C): 2 vials of recombinant Human IGSF4C		
	Assay Diluent A (Item D): 30 ml diluent buffer, 0.09% sodium azide as preservative. For		
	Standard/Sample (serum/plasma) diluent		
	<ul> <li>Assay Diluent B (Item E): 15 ml of 5x concentrated buffer. For Standard/Sample (cell culture medium/urine) diluent</li> </ul>		
	<ul> <li>Detection Affinity Reagent for IGSF4C (Item F): 2 vials of a 4x concentrated solution of anti- Human IGSF4C affinity reagent</li> </ul>		
	<ul> <li>IQELISA Detection Reagent (Item G): 1.4ml of a 10x concentrated stock</li> </ul>		
	Primer Solution (Item I): 1.7 ml vial		
	PCR Master Mix (Item J): 1.2 ml vial		
	PCR Preparation buffer (Item K): 1ml vial of 10x concentrated buffer		
	Final Wash Buffer (Item L): 10ml vial of 10x concentrated buffer		
Material not included:	Real-time PCR instrument, Bio-Rad recommended		
	<ul> <li>Precision pipettes to deliver 2 µL to 1 mL volumes</li> </ul>		
	Adjustable 1-25 mL pipettes for reagent preparation		
	<ul> <li>100 mL and 1 liter graduated cylinders</li> </ul>		
	Absorbent paper		
	Distilled or deionized water		
	<ul> <li>Log-log graph paper or computer and software for data analysis</li> </ul>		
	Tubes to prepare standard or sample dilutions		

• Heating block or water bath capable of 80°C

## Target Details

Target:	Cell Adhesion Molecule 4 (CADM4)
Alternative Name:	IGSF4C (CADM4 Products)
Gene ID:	199731
UniProt:	Q8NFZ8

## Application Details

Application Notes:	The Immuno-Quantitative ELISA (IQELISA) kits are an innovative assay platform that combines
	the specificity and ease of use of an ELISA with the sensitivity of real-time PCR. This results in
	an assay that is simultaneously familiar and cutting edge and enables the use of only 1/10th
	the sample volume while also providing 10x more sensitivity than a traditional ELISA.
Comment:	The Immuno-Quantitative ELISA (IQELISA™) kits are an innovative assay platform that
	combines the specificity and ease of use of an ELISA with the sensitivity of real-time PCR. Also

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN7526137 | 01/06/2025 | Copyright antibodies-online. All rights reserved.

	called immuno-PCR, this detection platform results in an assay that is simultaneously familiar and cutting edge. Compared to traditional ELISA, IQELISA™ enables the use of only 1/10th the
Sample Volume:	sample volume while also providing 10x more sensitivity. 25 μL
Plate:	Pre-coated
Protocol:	<ol> <li>Prepare all reagents, samples and standards as instructed</li> <li>Add 25 µL standard or sample to each well. Incubate for 2.5 hours at room temperature or overnight at 4°C</li> <li>Add 25 µL detection affinity reagent to each well. Incubate 1 hour at room temperature</li> </ol>
	4. Add 25µL of IQELISA Detection Reagent to each well. Incubate 1 hour at 100m temperature
	5. Add 15µL Primer solution 10µL of PCR master mix to each well
	6. Run real-time PCR
Restrictions:	For Research Use only
Handling	
Storage:	4 °C,-20 °C,-80 °C
Storage Comment:	May be stored for up to 6 months at 2° to 8°C from the date of shipment. Standard
	(recombinant protein) should be stored at -20°C or -80°C (recommended at -80°C) after
	reconstitution. Opened PCR plate or reagents may be stored for up to 1 month at 2° to 8°C.
	Note: the kit can be used within one year if the whole kit is stored at -20°C. Avoid repeated
	freeze-thaw cycles.
Expiry Date:	6 months