# antibodies .- online.com







### **OR1F1 Protein**



( )	ve	K\ /		A .
	$\cup$	1 V/	Щ.	V۷

Quantity:	50 μg
Target:	OR1F1
Origin:	Human
Source:	Mammalian Cells
Protein Type:	Synthetic Nanodisc

#### **Product Details**

Purpose:	Human OR1F1 full length protein-synthetic nanodisc	
Characteristics:	Unlike other membrane scaffold protein (MSP) Nanodisc on the market, our synthetic Nanodisc	
	can be prepared directly from the cells. The polymers used during this process have a dual	
	function. It dissolves the cell membranes, like the detergent, and uses cellular phospholipids to	
	form Nanodisc around the membrane proteins. The target protein embedded Nanodiscs can	
	then be purified.	

### **Target Details**

Target:	OR1F1
Alternative Name:	OR1F1 (OR1F1 Products)
Background:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated

## **Target Details**

	transduction of odorant signals. The olfactory receptor gene family is the largest in the	
	genome. The nomenclature assigned to the olfactory receptor genes and proteins for this	
	organism is independent of other organisms. [provided by RefSeq, Jul 2008]	
Molecular Weight:	The human full length OR1F1 protein has a MW of 34.9kDa	
UniProt:	043749	

Application Detai	ls
Comment:	Advantages of Synthetic Nanodiscs:
	Highly purified membrane proteins
	High solubility in aqueous solutions
	High stability
	<ul> <li>Proteins are in a native membrane environment and remain biologically active</li> </ul>
	<ul> <li>No detergent and can be used for cell-based assays</li> </ul>
	No MSP backbone proteins
	Limitations of Synthetic Nanodiscs:
	Intolerant to acids and high concentrations of divalent metal ions
Restrictions:	For Research Use only

## Handling

Format:	Lyophilized
Buffer:	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0).  Normally 5 % - 8 % trehalose is added as protectants before lyophilization.
Storage:	-20 °C,-80 °C
Storage Comment:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).  Lyophilized proteins are shipped at ambient temperature.
Expiry Date:	12 months