antibodies -online.com





Datasheet for ABIN7538435

OR2B3 Protein



Overview

Quantity:	50 μg
Target:	OR2B3
Origin:	Human
Source:	Mammalian Cells
Protein Type:	Synthetic Nanodisc

Product Details

Purpose:	Human OR2B3 full length protein-synthetic nanodisc
Characteristics:	Unlike other membrane scaffold protein (MSP) Nanodisc on the market, our synthetic Nanodisc
	can be prepared directly from the cells. The polymers used during this process have a dual
	function. It dissolves the cell membranes, like the detergent, and uses cellular phospholipids to
	form Nanodisc around the membrane proteins. The target protein embedded Nanodiscs can
	then be purified.

Target Details

Target:	OR2B3
Alternative Name:	OR2B3 (OR2B3 Products)
Background:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated

Target Details

	transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms.
Molecular Weight:	The human full length OR2B3 protein has a MW of 35.5 kDa
UniProt:	076000

Comment: Advantages of Synthetic Nanodiscs:	Application Details	
 Highly purified membrane proteins High solubility in aqueous solutions High stability Proteins are in a native membrane environment and remain biologically active No detergent and can be used for cell-based assays No MSP backbone proteins Limitations of Synthetic Nanodiscs: Intolerant to acids and high concentrations of divalent metal ions 	Comment:	 Highly purified membrane proteins High solubility in aqueous solutions High stability Proteins are in a native membrane environment and remain biologically active No detergent and can be used for cell-based assays No MSP backbone proteins Limitations of Synthetic Nanodiscs:

Restrictions: Fo

For Research Use only

Handling

Format:	Lyophilized
Buffer:	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5 % - 8 % trehalose is added as protectants before lyophilization.
Storage:	-20 °C,-80 °C
Storage Comment:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Expiry Date:	12 months