antibodies -online.com





Co to Droduct pag

Datasheet for ABIN7538438

OR2H2 Protein

Overview

| Quantity: | 50 µg |
|---------------|--------------------|
| Target: | OR2H2 |
| Origin: | Human |
| Source: | Mammalian Cells |
| Protein Type: | Synthetic Nanodisc |

Product Details

| Purpose: | Human OR2H2 full length protein-synthetic nanodisc |
|------------------|---|
| Characteristics: | Unlike other membrane scaffold protein (MSP) Nanodisc on the market, our synthetic Nanodisc |
| | can be prepared directly from the cells. The polymers used during this process have a dual |
| | function. It dissolves the cell membranes, like the detergent, and uses cellular phospholipids to |
| | form Nanodisc around the membrane proteins. The target protein embedded Nanodiscs can |
| | then be purified. |

Target Details

| Target: | OR2H2 |
|-------------------|---|
| Alternative Name: | OR2H2 (OR2H2 Products) |
| Background: | Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated |

Target Details

| | transduction of odorant signals. The olfactory receptor gene family is the largest in the |
|-------------------|---|
| | genome. The nomenclature assigned to the olfactory receptor genes and proteins for this |
| | organism is independent of other organisms. [provided by RefSeq, Jul 2008] |
| Molecular Weight: | The human full length OR2H2 protein has a MW of 34.8kDa |
| UniProt: | 095918 |

| UNIPROL: | 095918 |
|---------------------|--|
| Application Details | |
| Comment: | Advantages of Synthetic Nanodiscs: • Highly purified membrane proteins |
| | High solubility in aqueous solutions High stability |
| | Proteins are in a native membrane environment and remain biologically active No detergent and can be used for cell-based assays No MSP backbone proteins |
| | Limitations of Synthetic Nanodiscs: • Intolerant to acids and high concentrations of divalent metal ions |
| Restrictions: | For Research Use only |

Handling

| Format: | Lyophilized |
|------------------|--|
| Buffer: | Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5 % - 8 % trehalose is added as protectants before lyophilization. |
| Storage: | -20 °C,-80 °C |
| Storage Comment: | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Expiry Date: | 12 months |