antibodies -online.com





Datasheet for ABIN7538440

OR2M4 Protein



Overview

Quantity:	50 μg
Target:	OR2M4
Origin:	Human
Source:	Mammalian Cells
Protein Type:	Synthetic Nanodisc

Product Details

Purpose:

- 1	5 1 1 1 1 1 3 1 1 1 1 1 1 1 1 1 1 1 1 1
Characteristics:	Unlike other membrane scaffold protein (MSP) Nanodisc on the market, our synthetic Nanodisc
	can be prepared directly from the cells. The polymers used during this process have a dual
	function. It dissolves the cell membranes, like the detergent, and uses cellular phospholipids to
	form Nanodisc around the membrane proteins. The target protein embedded Nanodiscs can
	then be purified.

Human OR2M4 full length protein-synthetic nanodisc

Target Details

Target:	OR2M4
Alternative Name:	OR2M4 (OR2M4 Products)
Background:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response
	that triggers the perception of a smell. The olfactory receptor proteins are members of a large
	family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory
	receptors share a 7-transmembrane domain structure with many neurotransmitter and
	hormone receptors and are responsible for the recognition and G protein-mediated

Target Details

	transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]
Molecular Weight:	The human full length OR2M4 protein has a MW of 35.1kDa
UniProt:	Q96R27

Application Details

Com	me	nt:
-----	----	-----

Advantages of Synthetic Nanodiscs:

- · Highly purified membrane proteins
- · High solubility in aqueous solutions
- · High stability
- · Proteins are in a native membrane environment and remain biologically active
- · No detergent and can be used for cell-based assays
- No MSP backbone proteins

Limitations of Synthetic Nanodiscs:

· Intolerant to acids and high concentrations of divalent metal ions

Restrictions:

For Research Use only

Handling

Format:	Lyophilized
Buffer:	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5 % - 8 % trehalose is added as protectants before lyophilization.
Storage:	-20 °C,-80 °C
Storage Comment:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Expiry Date:	12 months