.-online.com antibodies

Datasheet for ABIN7538462 OR10H2 Protein



Overview

Quantity:	50 µg
Target:	OR10H2
Origin:	Human
Source:	Mammalian Cells
Protein Type:	Synthetic Nanodisc

Product Details

Purpose:	Human OR10H2 full length protein-synthetic nanodisc
Characteristics:	Unlike other membrane scaffold protein (MSP) Nanodisc on the market, our synthetic Nanodisc
	can be prepared directly from the cells. The polymers used during this process have a dual
	function. It dissolves the cell membranes, like the detergent, and uses cellular phospholipids to
	form Nanodisc around the membrane proteins. The target protein embedded Nanodiscs can
	then be purified.

Target Details

Target:	OR10H2
Alternative Name:	OR10H2 (OR10H2 Products)
Background:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response
	that triggers the perception of a smell. The olfactory receptor proteins are members of a large
	family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory
	receptors share a 7-transmembrane domain structure with many neurotransmitter and
	hormone receptors and are responsible for the recognition and G protein-mediated

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/2 | Product datasheet for ABIN7538462 | 06/14/2024 | Copyright antibodies-online. All rights reserved.

Target Details	
	transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]
Molecular Weight:	The human full length OR10H2 protein has a MW of 34.7kDa
UniProt:	060403

Application Details

Comment:	Advantages of Synthetic Nanodiscs:
	Highly purified membrane proteins
	High solubility in aqueous solutions
	High stability
	Proteins are in a native membrane environment and remain biologically active
	No detergent and can be used for cell-based assays
	No MSP backbone proteins
	Limitations of Synthetic Nanodiscs:
	Intolerant to acids and high concentrations of divalent metal ions
Restrictions:	For Research Use only
Handling	
Format:	Lyophilized
Buffer:	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0).
	Normally 5 % - 8 % trehalose is added as protectants before lyophilization.
Storage:	-20 °C,-80 °C
Storage Comment:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for
	use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).
	Lyophilized proteins are shipped at ambient temperature.
Expiry Date:	12 months