
# antibodies .- online.com







## **P2RY11 Protein**



( )	11/0	K\ /	iew	1
	$\cup$	ועוי	$\square \vee \vee$	ı

Quantity:	50 μg
Target:	P2RY11
Origin:	Human
Source:	Mammalian Cells
Protein Type:	Synthetic Nanodisc

#### **Product Details**

Purpose:	Human P2Y11 full length protein-synthetic nanodisc
Characteristics:	Unlike other membrane scaffold protein (MSP) Nanodisc on the market, our synthetic Nanodisc
	can be prepared directly from the cells. The polymers used during this process have a dual
	function. It dissolves the cell membranes, like the detergent, and uses cellular phospholipids to
	form Nanodisc around the membrane proteins. The target protein embedded Nanodiscs can
	then be purified.

### **Target Details**

Target:	P2RY11
Alternative Name:	P2Y11 (P2RY11 Products)
Background:	The product of this gene belongs to the family of G-protein coupled receptors. This family has several receptor subtypes with different pharmacological selectivity, which overlaps in some cases, for various adenosine and uridine nucleotides. This receptor is coupled to the stimulation of the phosphoinositide and adenylyl cyclase pathways and behaves as a selective
	purinoceptor. Naturally occuring read-through transcripts, resulting from intergenic splicing

### **Target Details**

	between this gene and an immediately upstream gene (PPAN, encoding peter pan homolog), have been found. The PPAN-P2RY11 read-through mRNA is ubiquitously expressed and encodes a fusion protein that shares identity with each individual gene product. [provided by RefSeq, Jul 2008]	
Molecular Weight:	The human full length P2Y11 protein has a MW of 40.3kDa	
UniProt:	Q96G91	
Pathways:	cAMP Metabolic Process	

#### **Application Details**

Co	-	5	-	~ +	٠.

Advantages of Synthetic Nanodiscs:

- · Highly purified membrane proteins
- · High solubility in aqueous solutions
- High stability
- · Proteins are in a native membrane environment and remain biologically active
- · No detergent and can be used for cell-based assays
- No MSP backbone proteins

Limitations of Synthetic Nanodiscs:

· Intolerant to acids and high concentrations of divalent metal ions

Restrictions:

For Research Use only

#### Handling

Format:	Lyophilized
Buffer:	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0).  Normally 5 % - 8 % trehalose is added as protectants before lyophilization.
Storage:	-20 °C,-80 °C
Storage Comment:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).  Lyophilized proteins are shipped at ambient temperature.
Expiry Date:	12 months