antibodies -online.com





Datasheet for ABIN7538586

VIPR1 Protein



Overview

Quantity:	50 μg
Target:	VIPR1
Origin:	Human
Source:	Mammalian Cells
Protein Type:	Synthetic Nanodisc

Product Details

Purpose:	Human VIPR1 full length protein-synthetic nanodisc
Characteristics:	Unlike other membrane scaffold protein (MSP) Nanodisc on the market, our synthetic Nanodisc
	can be prepared directly from the cells. The polymers used during this process have a dual
	function. It dissolves the cell membranes, like the detergent, and uses cellular phospholipids to
	form Nanodisc around the membrane proteins. The target protein embedded Nanodiscs can
	then be purified.

Target Details

Target:	VIPR1
Alternative Name:	VIPR1 (VIPR1 Products)
Background:	This gene encodes a receptor for vasoactive intestinal peptide, a small neuropeptide. Vasoactive intestinal peptide is involved in smooth muscle relaxation, exocrine and endocrine secretion, and water and ion flux in lung and intestinal epithelia. Its actions are effected through integral membrane receptors associated with a guanine nucleotide binding protein which
	activates adenylate cyclase. Several transcript variants encoding different isoforms have been

Target Details

	found for this gene. [provided by RefSeq, Oct 2011]
Molecular Weight:	The human full length VIPR1 protein has a MW of 51.5kDa
UniProt:	P32241

Molecular Weight:	The human full length VIPR1 protein has a MW of 51.5kDa
UniProt:	P32241
Application Details	
Comment:	Advantages of Synthetic Nanodiscs:
	 Highly purified membrane proteins High solubility in aqueous solutions High stability Proteins are in a native membrane environment and remain biologically active No detergent and can be used for cell-based assays No MSP backbone proteins Limitations of Synthetic Nanodiscs: Intolerant to acids and high concentrations of divalent metal ions
Restrictions:	For Research Use only
Handling	
Format:	Lyophilized
Buffer:	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5 % - 8 % trehalose is added as protectants before lyophilization.
Storage:	-20 °C,-80 °C
Storage Comment:	Store at -20°C to -80°C for 12 months in Ivonhilized form. After reconstitution, if not intended

. omiac.	
Buffer:	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5 % - 8 % trehalose is added as protectants before lyophilization.
Storage:	-20 °C,-80 °C
Storage Comment:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Expiry Date: 12 months