

Datasheet for ABIN7605575

anti-MMP13 antibody



Go to Product page

_					
	W	0	rv	10	W

Quantity:	100 μL
Target:	MMP13
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Monoclonal
Conjugate:	This MMP13 antibody is un-conjugated
Application:	Western Blotting (WB)

Product Details

Purpose:	Anti-MMP13 Rabbit Monoclonal Antibody
Immunogen:	A synthesized peptide derived from human MMP13
Clone:	IDA-13
Isotype:	IgG
Characteristics:	Anti-MMP13 Rabbit Monoclonal Antibody (ABIN7605575). Tested in WB application. This antibody reacts with Human, Mouse.
Purification:	Affinity-chromatography

Target Details

Target:	MMP13
Alternative Name:	MMP13 (MMP13 Products)

Target Details

Background:	Synonyms: Collagenase 3,3.4.24,Matrix metalloproteinase-13,MMP-13,MMP13, Tissue Specificity: Detected in fetal cartilage and calvaria, in chondrocytes of hypertrophic cartilage in vertebrae and in the dorsal end of ribs undergoing ossification, as well as in osteoblasts and periosteal cells below the inner periosteal region of ossified ribs. Detected in chondrocytes from in joint cartilage that have been treated with TNF and IL1B, but not in untreated chondrocytes. Detected in T lymphocytes. Detected in breast carcinoma tissue.
Molecular Weight:	52 kDa
UniProt:	P45452

Orm Foc.	
Application Details	
Application Notes:	WB 1:500-1:2000
Restrictions:	For Research Use only
Handling	
Format:	Liquid
Reconstitution:	Restore with deionized water (or equivalent) for reconstitution volume of 1.0 mL

Reconstitution.	Restore with deionized water (or equivalent) for reconstitution volume of 1.0 mL
Concentration:	Lot specific
Buffer:	Rabbit IgG in phosphate buffered saline, pH 7.4 , 150 mM NaCl, 0.02 % sodium azide and 50 % glycerol, 0.4 - 0.5 mg/mL BSA.
Preservative:	Sodium azide
Precaution of Use:	This product contains Sodium azide: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only.

Storage:	4 °C,-20 °C
Storage Comment:	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one
	month. Avoid repeated freeze-thaw cycles.