

# Datasheet for ABIN7637901

# anti-DNASEX antibody



#### Overview

Quantity:	100 μL
Target:	DNASEX
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This DNASEX antibody is un-conjugated
Application:	Western Blotting (WB), Immunohistochemistry (IHC), Immunoprecipitation (IP), Immunocytochemistry (ICC)

### **Product Details**

Purpose:	Polyclonal Antibody to Deoxyribonuclease X (DNASEX)
Immunogen:	RPE855Hu01Recombinant Deoxyribonuclease X (DNASEX)
Isotype:	IgG
Specificity:	The antibody is a rabbit polyclonal antibody raised against DNASEX. It has been selected for its ability to recognize DNASEX in immunohistochemical staining and western blotting.
Purification:	Antigen-specific affinity chromatography followed by Protein A affinity chromatography
Target Details	

Target:	DNASEX
Alternative Name:	DNASEX (DNASEX Products)

### **Target Details**

Background:	DNASE1L1, DNAS1L1, XIB, DNase-X, DNL1L, Deoxyribonuclease I-Like 1, Muscle-specific DNase	
	I-like	
UniProt:	P49184	
Application Details		
Application Notes:	Western blotting: 0.2-2 μg/mL,1:250-2500 Immunohistochemistry: 5-20 μg/mL,1:25-100	
	Immunocytochemistry: 5-20 $\mu$ g/mL,1:25-100 Optimal working dilutions must be determined by end user.	
Comment:	The thermal stability is described by the loss rate. The loss rate was determined by accelerated	
	thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious	
	degradation and precipitation were observed. The loss rate is less than 5% within the expiration	
	date under appropriate storage condition.	
Restrictions:	For Research Use only	
Handling		
Format:	Liquid	
Concentration:	500 μg/mL	
Buffer:	PBS, pH 7.4, containing 0.01 % SKL, 1 mM DTT, 5 % Trehalose and Proclin300.	
Preservative:	Dithiothreitol (DTT), ProClin	
Precaution of Use:	This product contains ProClin and Dithiothreitol (DTT): POISONOUS AND HAZARDOUS	
	SUBSTANCES which should be handled by trained staff only.	
Storage:	4 °C,-20 °C	
Storage Comment:	Store at 4°C for frequent use. Stored at -20°C in a manual defrost freezer for two year without detectable loss of activity. Avoid repeated freeze-thaw cycles.	