

# Datasheet for ABIN7638147

## anti-DNAH11 antibody



#### Overview

Quantity:	100 μL
Target:	DNAH11
Reactivity:	Human
Host:	Mouse
Clonality:	Monoclonal
Conjugate:	This DNAH11 antibody is un-conjugated
Application:	Immunohistochemistry (IHC), Western Blotting (WB), Immunoprecipitation (IP), Immunocytochemistry (ICC)

### **Product Details**

Target:

Alternative Name:

Purpose:	Monoclonal Antibody to Dynein, Axonemal, Heavy Chain 11 (DNAH11)
Immunogen:	RPJ201Hu02Recombinant Dynein, Axonemal, Heavy Chain 11 (DNAH11)
Clone:	C6
Specificity:	The antibody is a mouse monoclonal antibody raised against DNAH11. It has been selected for its ability to recognize DNAH11 in immunohistochemical staining and western blotting.
Purification:	Protein A + Protein G affinity chromatography
Target Details	

DNAH11

DNAH11 (DNAH11 Products)

## **Target Details**

Background:	DNAHBL, DNAHC11, DNHBL, DPL11, Dynein, Heavy Chain Beta-Like, Axonemal beta dynein
	heavy chain 11, Ciliary dynein heavy chain 11
UniProt:	Q96DT5
Application Details	
Application Notes:	Western blotting: 0.01-2 μg/mL,Immunohistochemistry: 5-20 μg/mL,Immunocytochemistry: 5-
	20 μg/mL,Optimal working dilutions must be determined by end user.
Comment:	The thermal stability is described by the loss rate. The loss rate was determined by accelerated
	thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious
	degradation and precipitation were observed. The loss rate is less than 5% within the expiration
	date under appropriate storage condition.
Restrictions:	For Research Use only
Handling	
Format:	Liquid
Concentration:	1 mg/mL
Buffer:	0.01M PBS, pH 7.4, containing 0.05 % Proclin-300, 50 % glycerol.
Preservative:	ProClin
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE which should be
	handled by trained staff only.
Storage:	4 °C,-20 °C
Storage Comment:	Store at 4°C for frequent use. Stored at -20°C in a manual defrost freezer for two year without
	detectable loss of activity. Avoid repeated freeze-thaw cycles.