

### Datasheet for ABIN7639500

# anti-GRIA4 antibody



#### Overview

Quantity:	100 μL
Target:	GRIA4
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This GRIA4 antibody is un-conjugated
Application:	Western Blotting (WB), Immunohistochemistry (IHC), Immunocytochemistry (ICC), Immunoprecipitation (IP)

### **Product Details**

Target:

Alternative Name:

GRIA4

**GRIA4 (GRIA4 Products)** 

Purpose:	Polyclonal Antibody to Glutamate Receptor, Ionotropic, AMPA 4 (GRIA4)
Isotype:	IgG
Specificity:	The antibody is a rabbit polyclonal antibody raised against GRIA4. It has been selected for its ability to recognize GRIA4 in immunohistochemical staining and western blotting.
Cross-Reactivity:	Mouse
Purification:	Antigen-specific affinity chromatography followed by Protein A affinity chromatography
Target Details	

### Target Details

Background:	GLUR4, GLUR4C, GLURD, AMPA-selective glutamate receptor 4
UniProt:	P48058
Pathways:	PI3K-Akt Signaling

## **Application Details**

Application Notes:	Western blotting: 0.2-2 μg/mL,1:250-2500 Immunohistochemistry: 5-20 μg/mL,1:25-100 Immunocytochemistry: 5-20 μg/mL,1:25-100 Optimal working dilutions must be determined by end user.
Comment:	The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration
Restrictions:	date under appropriate storage condition.  For Research Use only

## Handling

Format:	Liquid
Concentration:	500 μg/mL
Buffer:	PBS, pH 7.4, containing 0.01 % SKL, 1 mM DTT, 5 % Trehalose and Proclin300.
Preservative:	Dithiothreitol (DTT), ProClin
Precaution of Use:	This product contains ProClin and Dithiothreitol (DTT): POISONOUS AND HAZARDOUS SUBSTANCES which should be handled by trained staff only.
Storage:	4 °C,-20 °C
Storage Comment:	Store at 4°C for frequent use. Stored at -20°C in a manual defrost freezer for two year without detectable loss of activity. Avoid repeated freeze-thaw cycles.