

Datasheet for ABIN7640454

anti-HSD17B14 antibody



Overview

Quantity:	100 μL
Target:	HSD17B14
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This HSD17B14 antibody is un-conjugated
Application:	Western Blotting (WB), Immunohistochemistry (IHC), Immunocytochemistry (ICC), Immunoprecipitation (IP)

Product Details

Purpose:	Polyclonal Antibody to 17-Beta-Hydroxysteroid Dehydrogenase Type 14 (HSD17b14)
Isotype:	IgG
Specificity:	The antibody is a rabbit polyclonal antibody raised against HSD17b14. It has been selected for its ability to recognize HSD17b14 in immunohistochemical staining and western blotting.
Purification:	Antigen-specific affinity chromatography followed by Protein A affinity chromatography

Target Details

Target:	HSD17B14
Alternative Name:	HSD17b14 (HSD17B14 Products)
Background:	DHRS10, RetSDR3, SDR47C1, Retinal Short-Chain Dehydrogenase/Reductase 3, Short Chain
	Dehydrogenase/Reductase Family 47C,Member 1

Target Details

UniProt:	Q9BPX1
Application Details	
Application Notes:	Western blotting: 0.2-2 μg/mL,1:250-2500 Immunohistochemistry: 5-20 μg/mL,1:25-100 Immunocytochemistry: 5-20 μg/mL,1:25-100 Optimal working dilutions must be determined by end user.
Comment:	The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.
Restrictions:	For Research Use only
Handling	
Format:	Liquid
Concentration:	500 μg/mL
Buffer:	PBS, pH 7.4, containing 0.02 % Sodium azide, 50 % glycerol.
Preservative:	Sodium azide
Precaution of Use:	This product contains Sodium azide: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only.
Storage:	4 °C,-20 °C
Storage Comment:	Store at 4°C for frequent use. Stored at -20°C in a manual defrost freezer for two year without detectable loss of activity. Avoid repeated freeze-thaw cycles.