

### Datasheet for ABIN7641712

# anti-Cytokeratin 18 antibody



#### Overview

Quantity:	100 μL
Target:	Cytokeratin 18 (KRT18)
Reactivity:	Rat
Host:	Mouse
Clonality:	Monoclonal
Conjugate:	This Cytokeratin 18 antibody is un-conjugated
Application:	Western Blotting (WB), Immunohistochemistry (IHC), Immunoprecipitation (IP), Immunocytochemistry (ICC)

### **Product Details**

Purpose:	Monoclonal Antibody to Cytokeratin 18 (CK18)
Specificity:	The antibody is a mouse monoclonal antibody raised against CK18. It has been selected for its ability to recognize CK18 in immunohistochemical staining and western blotting.
Purification:	Antigen-specific affinity chromatography followed by Protein A affinity chromatography

## Target Details

Target:	Cytokeratin 18 (KRT18)
Alternative Name:	CK18 (KRT18 Products)
Background:	KRT18, CYK18, K18, Keratin 18, Cell proliferation-inducing gene 46 protein, Keratin, type I cytoskeletal 18
UniProt:	Q5BJY9

## **Target Details**

Pathways:	Apoptosis, Caspase Cascade in Apoptosis
Application Details	
Application Notes:	Western blotting: 0.2-2 $\mu$ g/mL,1:500-5000 Immunohistochemistry: 5-20 $\mu$ g/mL,1:50-200 Immunocytochemistry: 5-20 $\mu$ g/mL,1:50-200 Optimal working dilutions must be determined by end user.
Comment:	The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.
Restrictions:	For Research Use only
Handling	
Format:	Liquid
Concentration:	1 mg/mL
Buffer:	PBS, pH 7.4, containing 0.02 % Sodium azide, 50 % glycerol.
Preservative:	Sodium azide
Precaution of Use:	This product contains Sodium azide: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only.
Storage:	4 °C,-20 °C
Storage Comment:	Store at 4°C for frequent use. Stored at -20°C in a manual defrost freezer for two year without detectable loss of activity. Avoid repeated freeze-thaw cycles.