

Datasheet for ABIN7645527 **anti-CHN1 antibody**



[Go to Product page](#)

Overview

Quantity:	100 µL
Target:	CHN1 (ARHGAP2)
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This CHN1 antibody is un-conjugated
Application:	Western Blotting (WB), Immunohistochemistry (IHC), Immunocytochemistry (ICC), Immunoprecipitation (IP)

Product Details

Purpose:	Polyclonal Antibody to Chimerin 1 (CHN1)
Immunogen:	RPE318Hu01Recombinant Chimerin 1 (CHN1)
Isotype:	IgG
Specificity:	The antibody is a rabbit polyclonal antibody raised against CHN1. It has been selected for its ability to recognize CHN1 in immunohistochemical staining and western blotting.
Purification:	Antigen-specific affinity chromatography followed by Protein A affinity chromatography

Target Details

Target:	CHN1 (ARHGAP2)
Alternative Name:	CHN1 (ARHGAP2 Products)

Target Details

Background: CHN, DURS2, RhoGAP2, ARHGAP2, NC, N-Chimerin, Chimaerin, A-chimaerin, Alpha-chimerin, GTPase-Activating Protein,Rho 2, Duane Retraction Syndrome 2

UniProt: [P15882](#)

Application Details

Application Notes: Western blotting: 0.1-2 µg/mL,Immunohistochemistry: 5-20 µg/mL,Immunocytochemistry: 5-20 µg/mL,Optimal working dilutions must be determined by end user.

Comment: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Restrictions: For Research Use only

Handling

Format: Liquid

Concentration: 0.5 mg/mL

Buffer: 0.01M PBS, pH 7.4, containing 0.05 % Proclin-300, 50 % glycerol.

Preservative: ProClin, Sodium azide

Precaution of Use: This product contains ProClin and Sodium azide: POISONOUS AND HAZARDOUS SUBSTANCES which should be handled by trained staff only.

Storage: 4 °C,-20 °C

Storage Comment: Store at 4°C for frequent use. Stored at -20°C in a manual defrost freezer for two year without detectable loss of activity. Avoid repeated freeze-thaw cycles.