

## Datasheet for ABIN7646217

## anti-SLAMF7 antibody



O V CI V I C V V	
Quantity:	100 μL
Target:	SLAMF7
Reactivity:	Human
Host:	Mouse
Clonality:	Monoclonal
Conjugate:	This SLAMF7 antibody is un-conjugated
Application:	Western Blotting (WB), Immunohistochemistry (IHC), Immunoprecipitation (IP), Immunocytochemistry (ICC)
Product Details	
Purpose:	Monoclonal Antibody to Signaling Lymphocytic Activation Molecule Family, Member 7 (SLAMF7)
Immunogen:	RPK384Hu02Recombinant Signaling Lymphocytic Activation Molecule Family, Member 7 (SLAMF7)
Clone:	C2
Specificity:	The antibody is a mouse monoclonal antibody raised against SLAMF7. It has been selected for its ability to recognize SLAMF7 in immunohistochemical staining and western blotting.
Purification:	Protein A + Protein G affinity chromatography
Target Details	
Target:	SLAMF7

## **Target Details**

3	
Alternative Name:	SLAMF7 (SLAMF7 Products)
Background:	CD319, CRACC, 19A, CS1, CD2-like receptor-activating cytotoxic cells, Membrane protein FOAP-
UniProt:	Q9NQ25
Application Details	
Application Notes:	Western blotting: 0.5-2 μg/mL,Immunohistochemistry: 5-20 μg/mL,Immunocytochemistry: 5-20 μg/mL,Optimal working dilutions must be determined by end user.
Comment:	The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.
Restrictions:	For Research Use only
Handling	
Format:	Liquid
Concentration:	1 mg/mL
Buffer:	0.01M PBS, pH 7.4, containing 0.05 % Proclin-300, 50 % glycerol.
Preservative:	ProClin
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only.
Storage:	4 °C,-20 °C
Storage Comment:	Store at 4°C for frequent use. Stored at -20°C in a manual defrost freezer for two year without detectable loss of activity. Avoid repeated freeze-thaw cycles.