# ANTIBODIES ONLINE

# Datasheet for ABIN858555 TNFSF15 ELISA Kit

1 Image

4 Publications



Overview

Quantity:	96 tests
Target:	TNFSF15
Reactivity:	Mouse
Method Type:	Sandwich ELISA
Detection Range:	62.5-4000 pg/mL
Minimum Detection Limit:	62.5 pg/mL
Application:	ELISA

#### Product Details

Purpose:	This assay employs the quantitative sandwich enzyme immunoassay technique.
Sample Type:	Serum, Plasma, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Mouse TNFSF15.
Cross-Reactivity (Details):	Limited by current skills and knowledge, it is impossible for us to complete the cross-reactivity detection between the target antigen and all analogues for other species. Therefore, cross reaction may still exist.
Sensitivity:	15.6 pg/mL
Components:	<ul> <li>Assay plate (12 × 8 coated Microwells)</li> <li>Standard (freeze dried)</li> </ul>

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/6 | Product datasheet for ABIN858555 | 07/26/2024 | Copyright antibodies-online. All rights reserved.

- Biotin-antibody (100 × concentrate)
- HRP-avidin (100 × concentrate)
- Biotin-antibody Diluent
- HRP-avidin Diluent
- Sample Diluent
- Wash Buffer (25 × concentrate)
- TMB Substrate
- Stop Solution
- Adhesive Strip (for 96 wells)
- Instruction manual

### Target Details

Target:	TNFSF15
Alternative Name:	Tumor Necrosis Factor-Like Ligand 1 (TL1) (TNFSF15 Products)
Background:	Synonyms: MGC129934, MGC129935, TL1, TL1A, VEGI, VEGI192A, TNF ligand-related molecule 1 TNF superfamily ligand TL1A vascular endothelial cell growth inhibitor vascular endothelial growth inhibitor-192A
UniProt:	Q5UBV8
Pathways:	Positive Regulation of Endopeptidase Activity, Autophagy

## Application Details

Application Notes:	• The supplier is only responsible for the kit itself, but not for the samples consumed during the
Application Notes.	
	assay. The user should calculate the possible amount of the samples used in the whole test.
	Please reserve sufficient samples in advance.
	<ul> <li>Samples to be used within 5 days may be stored at 2-8°C, otherwise samples must be stored</li> </ul>
	at -20°C ( $\leq$ 1 month) or -80°C ( $\leq$ 2 months) to avoid loss of bioactivity and contamination.
	<ul> <li>Grossly hemolyzed samples are not suitable for use in this assay.</li> </ul>
	If the samples are not indicated in the manual, a preliminary experiment to determine the
	validity of the kit is necessary.
	• Please predict the concentration before assaying. If values for these are not within the range
	of the standard curve, users must determine the optimal sample dilutions for their particular
	experiments.
	Tissue or cell extraction samples prepared by chemical lysis buffer may cause unexpected
	ELISA results due to the impacts of certain chemicals.
	Owing to the possibility of mismatching between antigens from another resource and
	antibodies used in this supplier's kits (e.g., antibody targets conformational epitope rather
	than linear epitope), some native or recombinant proteins from other manufacturers may not
	be recognized by this supplier's products.

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/6 | Product datasheet for ABIN858555 | 07/26/2024 | Copyright antibodies-online. All rights reserved.

	<ul> <li>Influenced by factors including cell viability, cell number and cell sampling time, samples from cell culture supernatant may not be recognized by the kit.</li> </ul>
	<ul> <li>Fresh samples without long time storage are recommended for the test. Otherwise, protein degradation and denaturalization may occur in those samples and finally lead to wrong results.</li> </ul>
Comment:	Detection wavelength: 450 nm
	Information on standard material:
	Depending on the antigen to be detected, standards can be either native or recombinant
	protein. The recombinant proteins are being expressed in CHO cells in most cases. Please
	inquire for more information. The formulation of auxiliary material in the standard is considered
	proprietary information, however it does not contain any poisonous substance. Proclin 300
	(1:3000) is used as preservative.
	Information on reagents:
	In most cases the stop solution provided is 1 N H2SO4. The formulation of wash solution is
	proprietary information. None of the components contain (sodium) azide, thimerosal, 2-
	mercaptoethanol (2-ME) or any other poisonous materials. For the sandwich method kits, the
	sample diluent, antibody diluent, enzyme diluent and standard all contain BSA.
	Information on antibodies:
	The antibodies provided in different kits vary in regards to clonality and host. Some antibodies
	are affinity purified, some are Protein A
Sample Volume:	100 µL
Assay Time:	1 - 4.5 h
Plate:	Pre-coated
Protocol:	Antibody specific for TNFSF15 has been pre-coated onto a microplate. Standards and samples
	are pipetted into the wells and any TNFSF15 present is bound by the immobilized antibody.
	After removing any unbound substances, a biotin-conjugated antibody specific for TNFSF15 is
	added to the wells. After washing, avidin conjugated Horseradish Peroxidase (HRP) is added to
	the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate
	solution is added to the wells and color develops in proportion to the amount of TNFSF15


Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 3/6 | Product datasheet for ABIN858555 | 07/26/2024 | Copyright antibodies-online. All rights reserved.

measured.

<ul> <li>suck more than 10µL when pipetting.</li> <li>It is recommended to use distilled water to prepare reagents and samples. Using contaminated water or container for reagent preparation will influence detection result.</li> </ul>
It is recommended to use distilled water to prepare reagents and samples. Using
<ul> <li>Kindly use graduated containers to prepare the reagent. Please don't prepare the reagent directly in the Diluent vials provided in the kit.</li> <li>Bring all reagents to room temperature (18-25°C) before use for 30 min.</li> <li>Prepare fresh standard for each assay. Use within 4 hours and discard after use.</li> <li>Making serial dilution in the wells directly is not permitted.</li> <li>Please carefully reconstitute Standards according to the instruction. Avoid foaming and mix gently until the crystals have completely dissolved. To minimize imprecision caused by pipetting, use small volumes and ensure that pipettors are calibrated. It is recommended to avala mere then 10vl when pipetting.</li> </ul>
Note:
<ul> <li>HRP-avidin (1×) - Centrifuge the vial before opening. HRP-avidin requires a 100-fold dilution. The suggested dilution is 10µL of HRP-avidin + 990µL of HRP-avidin Diluent.</li> <li>Wash Buffer (1×) - If crystals have formed in the concentrate, warm up to room temperature and mix gently until the crystals have completely dissolved. Dilute 20mL of Wash Buffer (1×).</li> <li>Standard - Centrifuge the standard vial at 6000-10000rpm for 30s. Reconstitute the Standard with 1ml of Sample Diluent. Do not substitute other diluents. This reconstitution produces a stock solution. Mix the standard to ensure complete reconstitution and allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.</li> <li>Pipette 250µL of Sample Diluent into each tube. Use the stock solution to produce a 2-fold dilution series. Mix each tube thoroughly before the next transfer. The undiluted Standard serves as the high standard. Sample Diluent serves as the zero standard (0ng/mL).</li> </ul>

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 4/6 | Product datasheet for ABIN858555 | 07/26/2024 | Copyright antibodies-online. All rights reserved.

	protection when using this material.
Handling Advice:	<ul> <li>The kit should not be used beyond the expiration date on the kit label.</li> <li>Do not mix or substitute reagents with those from other lots or sources.</li> <li>If samples generate values higher than the highest standard, dilute the samples with Sample</li> </ul>
	Diluent and repeat the assay.
	<ul> <li>Any variation in Sample Diluent, operator, pipetting technique, washing technique, incubation time/temperature and kit age can cause variation in binding.</li> </ul>
	This assay is designed to eliminate interference by soluble receptors, binding proteins and
	other factors present in biological samples. Until all factors have been tested in the
	Immunoassay, the possibility of interference cannot be excluded.
Storage:	4 °C/-20 °C
Storage Comment:	For unopened kit: All the reagents should be kept according to the labels on vials.
Expiry Date:	6 months
Publications	
Product cited in:	Pham, ODonnell, Al-Shamkhani, Kerrinnes, Tsolis, McSorley: "T cell expression of IL-18R and
	DR3 is essential for non-cognate stimulation of Th1 cells and optimal clearance of intracellular
	bacteria." in: PLoS pathogens, Vol. 13, Issue 8, pp. e1006566, (2017) (PubMed).
	Qin, Xu, Qi, Yang, Zhang, Liu, Xu, Xiang, Xiao, Cao, Wei, Zhang, Li: "Tumour necrosis factor
	superfamily member 15 (Tnfsf15) facilitates lymphangiogenesis via up-regulation of Vegfr3
	gene expression in lymphatic endothelial cells." in: The Journal of pathology, Vol. 237, Issue 3,
	pp. 307-18, (2015) (PubMed).
	Liu, Wang, Ming, Song, Zhang, Chen, Zeng, Mao: "S100A9: A Potential Biomarker for the
	Progression of Non-Alcoholic Fatty Liver Disease and the Diagnosis of Non-Alcoholic
	Steatohepatitis." in: <b>PLoS ONE</b> , Vol. 10, Issue 5, pp. e0127352, (2015) (PubMed).
	Voloshin, Veselski?: "[Lipid peroxidation in patients with the initial stage of circulatory
	encephalopathy of atherosclerotic origin]." in: <b>Vrachebnoe delo</b> , Issue 9, pp. 90-1, (1990) (



#### ELISA

Image 1. Typical standard curve

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 6/6 | Product datasheet for ABIN858555 | 07/26/2024 | Copyright antibodies-online. All rights reserved.