# antibodies .- online.com





## anti-Amylin/DAP antibody (Alexa Fluor 647)



Go to Product page

$\sim$			
	N/P	r\/I	i⊢₩

Quantity:	100 μL	
Target:	Amylin/DAP (IAPP)	
Reactivity:	Human, Rat, Mouse	
Host:	Rabbit	
Clonality:	Polyclonal	
Conjugate:	This Amylin/DAP antibody is conjugated to Alexa Fluor 647	
Application:	Western Blotting (WB), Immunofluorescence (Paraffin-embedded Sections) (IF (p))	

#### **Product Details**

Immunogen:	KLH conjugated synthetic peptide derived from human Amylin	
Isotype:	IgG	
Cross-Reactivity:	Human, Mouse, Rat	
Purification:	Purified by Protein A.	

### Target Details

Target:	Amylin/DAP (IAPP)	
Alternative Name:	Amylin (IAPP Products)	
Background:	Synonyms: Islet amyloid polypeptide, IAPP, DAP, Diabetes associated peptide, IAP, Insulinoma amyloid peptide.	
	Background: Amylin is commonly found in pancreatic islets of patients suffering diabetes	
	mellitus type II, or harboring an insulinoma. While the assosciation of amylin with the	

#### **Target Details**

Expiry Date:

12 months

rarget Details		
	development of type II diabetes has been known for some time, a direct causative role for amylin has been harder to establish. Recent results suggest that amylin, like the related beta	
	amyloid (Abeta) assosciated with Alzheimer's disease, can induce apoptotic cell death in particular cultured cells, an effect that may be relevant to the development of type II diabetes.	
Gene ID:	3375	
Pathways:	Hormone Activity, Feeding Behaviour	
Application Details		
Application Notes:	IF(IHC-P): (1:50-200)	
	Optimal working dilution should be determined by the investigator.	
Restrictions:	For Research Use only	
Handling		
Format:	Liquid	
Concentration:	1 μg/μL	
Buffer:	Aqueous buffered solution containing 0.01M TBS ( pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.	
Preservative:	Sodium azide	
Precaution of Use:	This product contains Sodium azide: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only.	
Storage:	-20 °C	
Storage Comment:	Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.	