



[Go to Product page](#)

Datasheet for ABIN888020
anti-CCDC65 antibody (AA 301-400) (Cy7)

Overview

| | |
|----------------------|---------------------------------------------------------------------------------------------------------|
| Quantity: | 100 µL |
| Target: | CCDC65 |
| Binding Specificity: | AA 301-400 |
| Reactivity: | Human |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Conjugate: | This CCDC65 antibody is conjugated to Cy7 |
| Application: | Immunofluorescence (Cultured Cells) (IF (cc)), Immunofluorescence (Paraffin-embedded Sections) (IF (p)) |

Product Details

| | |
|-----------------------|------------------------------------------------------------|
| Immunogen: | KLH conjugated synthetic peptide derived from human CCDC65 |
| Isotype: | IgG |
| Predicted Reactivity: | Human,Mouse,Rat,Dog,Cow,Sheep,Pig,Horse,Rabbit |
| Purification: | Purified by Protein A. |

Target Details

| | |
|-------------------|--------------------------------------------------------------------------------------|
| Target: | CCDC65 |
| Alternative Name: | CCDC65 (CCDC65 Products) |
| Background: | Synonyms: CCD65_HUMAN, CCDC65, Coiled coil domain containing 65, Coiled-coil domain- |

Target Details

containing protein 65, NYD-SP28, Testis development protein NYD-SP28.

Background: CCDC65 is a sperm tail protein that is posttranslationally modified during sperm capacitation (Zheng et al., 2006 [PubMed 17089017]).[supplied by OMIM, Mar 2008].

Gene ID: 85478

Application Details

Application Notes: IF(IHC-P) 1:50-200
IF(IHC-F) 1:50-200
IF(ICC) 1:50-200

Restrictions: For Research Use only

Handling

Format: Liquid

Concentration: 1 µg/µL

Buffer: Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.

Preservative: ProClin

Precaution of Use: This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.

Storage: -20 °C

Storage Comment: Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.

Expiry Date: 12 months