antibodies -online.com





anti-CYP46A1 antibody (AA 51-150) (AbBy Fluor® 555)



Go to Product page

\sim					
()	VE	۲۱	/1	\triangle	Λ

Quantity:	100 μL	
Target:	CYP46A1	
Binding Specificity:	AA 51-150	
Reactivity:	Mouse	
Host:	Rabbit	
Clonality:	Polyclonal	
Conjugate:	This CYP46A1 antibody is conjugated to AbBy Fluor® 555	
Application:	Western Blotting (WB), Immunofluorescence (Cultured Cells) (IF (cc)), Immunofluorescence (Paraffin-embedded Sections) (IF (p))	

Product Details

Immunogen:	KLH conjugated synthetic peptide derived from human CYP46	
Isotype:	IgG	
Cross-Reactivity:	Mouse	
Predicted Reactivity:	y: Human,Rat,Dog,Cow,Sheep,Horse,Rabbit	
Purification:	Purified by Protein A.	

Target Details

Target:	CYP46A1
Alternative Name:	CYP46 (CYP46A1 Products)

Target Details

Background:	Synonyms: Cholesterol 24-hydroxylase, CP46, CYP46, CP46A_HUMAN, cytochrome P450 family		
	46.		
	Background: The Cyp46 gene encodes a member of the cytochrome P450 superfamily of		
	enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions		
	involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This		
	endoplasmic reticulum protein is expressed exclusively in the brain, where it converts		
	cholesterol to 24S-hydroxycholesterol by adding a hydroxyl group to cholesterol, producing a		
	product that is more soluble than cholesterol and able to be exported from the brain.Cyp46 is		
	also known as 24S-cholesterol hydroxylase.		
Gene ID:	10858		
Application Details			
Application Notes:	IF(IHC-P) 1:50-200		
	IF(IHC-F) 1:50-200		
	IF(ICC) 1:50-200		
Restrictions:	For Research Use only		
Handling			
Format:	Liquid		
Concentration:	1 μg/μL		
Buffer:	Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and		
	50 % Glycerol.		
Preservative:	ProClin		
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be		
	handled by trained staff only.		
Storage:	-20 °C		
Storage Comment:	Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.		
Expiry Date:	12 months		