

[Go to Product page](#)

Datasheet for ABIN903928  
**anti-KLHL20 antibody (AA 13-100) (Cy5)**

## Overview

Quantity:	100 µL
Target:	KLHL20
Binding Specificity:	AA 13-100
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This KLHL20 antibody is conjugated to Cy5
Application:	Immunofluorescence (Cultured Cells) (IF (cc)), Immunofluorescence (Paraffin-embedded Sections) (IF (p))

## Product Details

Immunogen:	KLH conjugated synthetic peptide derived from human KLHL20
Isotype:	IgG
Cross-Reactivity:	Rat
Predicted Reactivity:	Human, Mouse, Cow, Sheep, Chicken, Rabbit
Purification:	Purified by Protein A.

## Target Details

Target:	KLHL20
Alternative Name:	KLHL20 ( <a href="#">KLHL20 Products</a> )

## Target Details

Background:	<p>Synonyms: Kelch like ECT2 interacting protein, Kelch like protein 20, Kelch like protein X, Kelch mot containing protein, Kelch-like ECT2-interacting protein, Kelch-like protein 20, Kelch-like protein X, KHLHX, KLEIP, KLH20CHUMAN, KLHL20, KLHLX, RP3-383J4.3.</p> <p>Background: Substrate-specific adapter of a BCR (BTB-CUL3-RBX1) E3 ubiquitin-protein ligase complex involved in interferon response. The BCR(KLHL20) E3 ubiquitin ligase complex mediates the ubiquitination of DAPK1, leading to its degradation by the proteasome, thereby acting as a negative regulator of apoptosis. Also acts as a regulator of endothelial migration during angiogenesis by controlling the activation of Rho GTPases.</p>
Gene ID:	27252

## Application Details

Application Notes:	<p>IF(IHC-P) 1:50-200</p> <p>IF(IHC-F) 1:50-200</p> <p>IF(ICC) 1:50-200</p>
Restrictions:	For Research Use only

## Handling

Format:	Liquid
Concentration:	1 µg/µL
Buffer:	Aqueous buffered solution containing 0.01M TBS ( pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.
Preservative:	ProClin
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.
Storage:	-20 °C
Storage Comment:	Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.
Expiry Date:	12 months