

Datasheet for ABIN955943

Peroxidase Labeling -NH₂ Kit



[Go to Product page](#)

Overview

Quantity: 3 tests

Target: Peroxidase Labeling-NH₂

Application: Labeling (LbI)

Product Details

Components: Peroxidase, NH₂-reactive 3 X 100 µg
Reaction Buffer 200 µL
Wash Buffer 4 mL
.

Material not included: 0.5 mL microtubes.
10 µL and 200 µL adjustable pipettes.
Microcentrifuge
37°C Incubator

Target Details

Target: Peroxidase Labeling-NH₂

Application Details

Restrictions: For Research Use only

Handling

Storage: 4 °C

Handling

Storage Comment:	Store all components at 4°C. Stable for 6 months at 4°C with protection from moisture. Storage Buffer 4 mL
Expiry Date:	6 months
