

Datasheet for ABIN956023

Peroxidase Labeling - SH Kit



[Go to Product page](#)

Overview

Quantity:	3 tests
Target:	Peroxidase Labeling-SH
Application:	Labeling (LbI)

Product Details

Components:	Peroxidase, SH-reactive 3 X 100 µg Reducing Agent 3 tubes Solution A 4 mL Solution B 1 mL Reaction Buffer 200 µL .
Material not included:	0.5 mL microtubes. 10 µL and 200 µL adjustable pipettes. Microcentrifuge 37°C Incubator

Target Details

Target:	Peroxidase Labeling-SH
---------	------------------------

Application Details

Restrictions:	For Research Use only
---------------	-----------------------

Handling

Storage:	4 °C
Storage Comment:	Store all components at 4°C. Stable for a year at 4°C with protection from moisture. Storage Buffer 4 mL
Expiry Date:	12 months