.-online.com antibodies

Datasheet for ABIN958833 anti-Nuclease S1 antibody (Biotin)



Overview	
Quantity:	1 mL
Target:	Nuclease S1 (NUCS)
Reactivity:	Aspergillus oryzae
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This Nuclease S1 antibody is conjugated to Biotin
Application:	Immunofluorescence (IF), Immunohistochemistry (IHC), ELISA, Western Blotting (WB)
Product Details	
Immunogen:	Aspergillus oryzae Nuclease S1
	Type of Immunogen: Purified protein
Specificity:	Aspergillus oryzae Nuclease S1. The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrop horesis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and enzyme inhibition. Cross-reactivity Cross-reactivities against enzymes of other sources may occur but have not been determined.
Purification:	Ion exchange chromatography
Target Details	
Target:	Nuclease S1 (NUCS)

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/2 | Product datasheet for ABIN958833 | 09/12/2023 | Copyright antibodies-online. All rights reserved.

Target Details	
Alternative Name:	Nuclease S1 (NUCS Products)
Application Details	
Application Notes:	Approved: ELISA, IF, IHC, WB
	Usage: This product is intended for use in precipitating and non-precipitating antibody-binding assays (such as e.g., ELISA and Western blotting and immunofluorescence or histochemical techniques). The applications listed have been tested for the unconjugated form of this product. Other forms have not been tested.
Restrictions:	For Research Use only
Handling	
Format:	Lyophilized
Reconstitution:	Sterile distilled water 1 ml
Concentration:	Lot specific
Buffer:	Lyophilized from PBS, pH 7.2
Handling Advice:	Avoid repeat freeze-thaw cycles.
Storage:	4 °C,-20 °C
Storage Comment:	Long term: -20°C Short term: +4°C. Avoid repeat freeze-thaw cycles.