

Datasheet for ABIN965411

Western Blot Kit Chemiluminescent Western Blot Kit for use with Rabbit Primary Antibody FemtoMax™



[Go to Product page](#)

2 Images

Overview

Quantity:	1 kit
Application:	Western Blotting (WB)

Product Details

Brand:	FemtoMax™
--------	-----------

Target Details

Background:	<p>This Chemiluminescent Western Blot Kit allows for the detection of primary Rabbit polyclonal or monoclonal antibodies provided by the user. After protein separation and transfer, the membrane is probed with primary antibody. Detection of the membrane bound primary antibody-antigen complex is achieved by the addition of a secondary antibody conjugated to the enzyme horseradish peroxidase. The enzyme reacts with a specialized formulation of luminol, an extremely sensitive, non-radioactive substrate that emits light and allows visualization using X-ray film or other imaging methods, including highly sensitive CCD cameras and imaging systems. Because of the extremely high sensitivity of our FemtoMax[®] Chemiluminescent Substrate Kit for Western Blotting, primary and secondary antibodies can be used at greater dilutions. If you currently dilute your primary antibody 1:1,000 using an ECL substrate equivalent, then dilute your antibody 1:5,000 using the FemtoMax[®] substrate.</p> <p>Synonyms: Western Blotting Kit, Chemiluminescent Kit, Peroxidase Kit, Immunoblotting Kit</p>
-------------	--

Application Details

Application Notes:	<p>Chemiluminescent Western Blot Kit, namely our FemtoMax™ Chemiluminescent Substrate Kit for Western Blotting, combines all of the necessary reagents with a rapid proven protocol and the extremely high signal detection of our luminol substrate. The Chemiluminescent Western</p>
--------------------	--

Application Details

Blot Kit design includes straightforward procedures and color-coding to allow for ease of use. This kit contains sufficient substrate for up to 30 mini blots at 7.5 x 8 cm² and is stable for at least 1 year when stored at +4°C.

Comment: Detection Kit Type: Chemiluminescent Western Blot Kit


Restrictions: For Research Use only

Handling

Handling Advice: Wash buffers MUST NOT contain SODIUM AZIDE or other inhibitors of peroxidase activity!

Storage: RT/4 °C

Images



Western Blotting

Image 1.



Western Blotting

Image 2. FemtoMax™ Chemiluminescent Western Blotting kit delivers: Super Sensitive Signal. Conserves Costly Antibody. Femtogram Sensitivity. Excellent Stability. Platform Flexibility.